

The C2 CryoBalloon™ Ablation System

For the Endoscopic Cryoablation of Barrett's Oesophagus
and Esophageal Squamous Dysplasia




Deep Controlled Ablation in a Single Stroke

The NEW C2 Cryoballoon Ablation System combines predictable dosing with ablation extending to submucosa that only Cryoballoon can provide.



Ablation extending to the submucosa¹



Combines deep ablation with a potentially favorable safety profile²

Effective Endoscopic Treatment For BE-Associated Neoplasia²

92%

Complete Eradication of Dysplasia

84%

Complete Eradication of Intestinal Metaplasia

Less Pain⁴

PATIENTS REPORTING SEVERE PAIN AT DAY 2

RFA

CryoBalloon

% 0 5 10 15 20 25 30 35 40 45 50

1. Louis B, et al. Evaluation of a Novel Swiggy Cryoballoon Ablation System in Bench, Porcine, and Human Esophageal Models. Gastrointestinal Endoscopy, Volume 85, No. 5S, AB513-AB514; 2017. 2. Schwartz D, et al. Endoscopic Cryotherapy of Barrett's Esophagus With a Novel Local Ablation Device: Feasibility and Safety in Porcine and Human Models. Techniques in Gastrointestinal Endoscopy, 2017 Nov (2011); 8. 3. Manca Irene Camis, Et Al. Multichannel Nitrous Oxide Cryoballoon Ablation With or Without Endoscopic Mucosal Resection (Emr) For Treatment of Neoplastic Barrett's Esophagus: Preliminary Results of a Prospective Clinical Trial in Treatment-Naïve and Previously Ablated Patients. P556, United European Gastroenterology Journal, 2016; 2 (Supplement 1). 4. Van Munster, et al. Post-Procedural Pain After Endoscopic Ablation Therapy of Barrett's Esophagus: Is Cryo Ablation Less Painful Compared with Radiofrequency Ablation?. Gastrointestinal Endoscopy, 2017 May, Volume 85, Issue 5, Supplement. Page AB088.



Pear Ablation Catheter
Designed for Ablation of GEJ



Standard Shaped Catheter
for Tubular Esophagus

All in a Smart, Economic + Compact Design

C2 CryoBalloon Ablation System Provides:

- Targeted access delivered through the working channel of a therapeutic endoscope
- An intuitive interface
- Smart, compact design that eliminates the need for large tanks or capital equipment



C2 CryoBalloon Catheters

- Multiple shapes & targeting options designed for various anatomic structures



C2 CryoBalloon Controller

- Easy to program dosing
- Houses C2 CryoBalloon cartridges



C2 CryoBalloon Foot Pedal

- Enables rotational & axial targeting

VARIMED

Dział Sprzedaży / Akcesoria:
Varimed Sp. z o.o.
Wrocław 50-442
ul. T. Kościuszki 115/4U
tel. +48 71 361 14 18
fax: +48 71 361 86 11

Dział Sprzedaży:
Leszno 64-100
ul. J. Karasia 16
tel. +48 65 529 36 89
fax: +48 65 520 99 78
e-mail: varimed@varimed.pl

Dział Serwisu:
Leszno 64-100
ul. Uslugowa 20
tel. +48 65 529 72 35
fax: +48 65 529 36 97
e-mail: serwis@varimed.pl